



SMARTCARE
Presented by WENETIS ENTERPRISES

Best Asset Management Software





Biomedical Assets

The software build with the experience of biomedical engineers' guidance and tested with actual data.

BME SmartCare follows World standards with WHO guidelines

We also customize for other departments



IT Assets



Facility Assets



SMARTCARE SOFTWARE



SMARTCARE
Presented by WENETIS ENTERPRISES

- ✓ Asset Management
- ✓ Service Call Management
- ✓ Cost Analysis
- ✓ Multi Location Hospital Management
- ✓ Spare Part Management
- ✓ Reports
- ✓ PM & Calibration Scheduling

- Graphical Dashboard ✓
- SERVICE CALL & PM GRAPH ✓
- Cost Expenses Percentage ✓
- One Click Equipment History ✓
- Assign Service Calls Facility ✓



DASHBOARD



REPORTING

- ✓ Equipment Details Report
- ✓ MIS Report
- ✓ Part Used Report
- ✓ Uptime, Downtime & Response Time Report
- ✓ MTR & MTBF Report
- ✓ Equipment Comparison Reports

- Barcode & Qr Code Scanning ✓
- Manage Preventive Maintenance ✓
- Service Call Management ✓
- Live Chat Service ✓
- Digital Signature On Report ✓
- 100% Paperless ✓



ANDROID APPLICATION



INFINITE STUFF



SMARTCARE

- ✓ Unlimited Users
- ✓ Unlimited Locations/Customer
- ✓ Unlimited Equipment's
- ✓ Unlimited Cloud Space
- ✓ Unlimited Records

- AMC/CMC Notification ✓
- PM/Calibration Notification ✓
- Warranty Expire Notification ✓
- License Notification ✓
- Insurance & QA Notification ✓
- Extra Cost Expenses Notification ✓
- New Service Call Notification ✓



NOTIFICATION



BENEFITS

- ✓ Structurize Data
- ✓ Auto PM/Calibration Scheduling
- ✓ Save Time
- ✓ Save Money
- ✓ On One Click Records Availability
- ✓ Easy To Use
- ✓ Engineer Work Analysis
- ✓ Scan To View Asset Details
- ✓ Inventory Management With Cost
- ✓ Auto Email & SMS To Vendor/Customer

- Multiple Department ✓
- Multiple Hospital/clinic ✓
- Multiple Location ✓
- Centralize Dashbaord ✓



MULTI DEPARTMENT & MULTI LOCATION


Best Asset Management Software




Screen Shots



Customizable Dashboard



List of Preventive Maintenance



List of Preventive Maintenance

Dashboard / Preventive Maintenance

Search Maintenance Date :




Search

+ Add new Preventive Maintenance

Export to Excel

Show 10 entries

Search:



Equipment Name	Maintenance DateTime	Company Engineer Name	Biomedical Engineer Name	Part Details & condition	Next Pm Due Date	Created Date	Is Enable	Created By	Action
1	11/26/2016 12:00:00 AM	Pratik Panbude	Haridas More	Test	11/26/2016 12:00:00 AM	11/26/2016 3:40:20 PM	Active	1	  

Showing 1 to 1 of 1 entries

Previous Next



Create PM list


 

Create

[Dashboard](#) / Create


<p>Equipment Name</p> <input type="text" value="Select--"/>	<p>Part Details & condition</p> <input type="text"/>
<p>Maintainance DateTime</p> <input type="text"/>	<p>Comments</p> <input type="text"/>
<p>Check List</p> <input type="button" value="Choose File"/> No file chosen	<p>Feedback Form</p> <input type="text"/>
<p>PM Number</p> <input type="text"/>	
<p>Company Engineer Name</p> <input type="text"/>	
<p>Biomedical Engineer Name</p> <input type="text"/>	

PM Report

2017 © Application by Wenetis. 

Service call creation

3 2 2

 admin

Select--

Call Time

Call Urgency

Low

Complaint Number

Description

Eye Witness

Head Of Department

Bio-Medical Engineer

Remark

Solved By

InHouse Company

Parts Used

Yes No

Report Copy

No file chosen

Rating

☆☆☆☆☆

Status

Select

2017 © Application by Wenetis.

Service call entry fields

Search By
 Asset Code By Name

Asset Code

Complaint By

Complaint Type

Critical / Non-Critical

Complaint Time

Call Attended Time

Description

Service Call Solution

Service Call Cost

Cost Saved

Upload Files
 No file chosen

Solved By
 InHouse Company

Remark

BioMedical Engineer

Status

Complaint Closed By

Complaint Close Time

Date	Followup	Action
------	----------	--------

Parts Used
 Yes No

Approval Needed
 Yes No

© smartcars@wonele.com | +91 9623813114 | EME SmartCars Developed by Wonele Enterprise

Service call report

Joblog Note

Complaint Number	Biomedical/22-12-2020/2	Job Date	22-12-2020 11:00
Complaint Type	AMC	Critical/Non-Critical	Critical
Job done by	Yogesh Mhase	Attended On	22-12-2020 11:05
Response Time	5	Job completed Time	22-12-2020 11:20
Problem	Image not display on monitor.		
Action Taken	Checked device and reinstalled connections. Now working ok.(But machine under observation) Checked by company. Multicore cable needs to be replace.(23-12-2020)		
Service call cost	0		
Remark			

Equipment Name	C-Arm	Asset Code	2115
Manufacturer	Surgico-60D-HF	Department Name	NEURO OT
Model	Phillips	Serial No	G8L211
Date Of Installation	24-11-2008	Contract Details	AMC
PO Number		Purchase Cost	740000
Equipment Expiry Date	22-11-2018	Equipment Used	122.52 %
Complaints	6	Equipment Life	12 Years 2 Months 28 Days
UpTime	99.94 %	DownTime	63.45 Hrs
Warranty Due Date	24-11-2009	Warranty In Years	1

Follow up Details

Date	Action Taken
23-12-2020	Checked by company .Problem found in the multicore cable .Video signal cable found internal open. Now temporary connected non shielded wire .Multicore cable needs to be replace.
10-02-2021	replaced Multicore cable

Used Part Details

Name	Quantity	Price
Multicor cable	1	17500.00

Bio-medical Engineer

Users Sign

Attached Files

- c-arm-g8l112.jpg
- 032-OHR-JAN-Deccan-Health-Care.pdf
- C-arm-multicore-cable-1-.jpg

Print

Cancel

End-User complaints entry

Add Breakdown

Asset Code

Action

Equipment

Critical / Non-Critical

Complaint Type

Complaint By

Problem

Upload Files(Image/Video/Files)

 No file chosen

PM report

Preventive Maintenance Details

(Basic information)

Equipment Name	Syringe pump	Model Name	SP104
Serial No	1310135	Manufacturer	Skandray
Asset Code	893	Maintenance Date	04-01-2020
Company Engineer Name	Test	Biomedical Engineer Name	Biomedical Admin
PartDetails And Condition		PM Scheduled Date	01-03-2019
Feedback Form		NextPM Due Date	04-07-2020
PM CheckList		PM Report	

Bio-medical Engineer

Users Sign





Equipment history-1

Select Name-Asset Code-Model No-Serial No- - -

Search

Equipment Details

Export all to excel

Equipment Details			
Equipment Name :	Ventilator	AssetCode :	1
Hospital :	Biomedical	Building :	Hightech Building
Floor :	Fifth Floor	Department :	MICU
Total Complaints :	10	Manufacturer :	Maquet
ResponceTime In Min:	9	SerialNo :	37861
Installation Date :	23-11-2016	Equipment Model :	servo"s
Equipment Expiry :	4 Years 2 Months 29 Days	DownTime :	1611
Purchase Cost :	761905	Cost Expenditure :	3.88 %
Actions	   		

Equipment history-2

Service Call Details																			
Action	Status	Complaint Number	AssetCode	Equipment Name	Model	Serial No	Manufacturer	Department Name	Problem	Action Taken	CallDateTime	CallClosedTime	DownTime (Hrs. Min)	ResponseTime InMin	Service Cost	Cost Saved	Biomedical Engineer	IsEnable	Location
	Call Closed	Biomedical/20-03-2020/1	1	Ventilator	servo's	37861	Maquet	MICU	Not working . power failure	Solvedd	20-03-2020 09:10	20-03-2020 01:00	0.50	25	1000.00	5000.00	Biomedical Admin	Active	Biomedical
	Call Closed	Biomedical/30-03-2020/1	1	Ventilator	servo's	37861	Maquet	MICU	internal leakage	Pre use check test done.	30-03-2020 14:29		0.00	1	0.00	0.00	Manisha mufe	Active	Biomedical
	Call Closed	Biomedical/09-06-2020/3	1	Ventilator	servo's	37861	Maquet	MICU	pre use check failed		09-06-2020 16:00	09-06-2020 16:15	0.15	10	0.00	0.00	Manisha mufe	Active	Biomedical
	Call Closed	Biomedical/09-06-2020/4	1	Ventilator	servo's	37861	Maquet	MICU	internal leakage	Pre use check test passed . machine working ok.	09-06-2020 16:28	09-06-2020 16:45	0.17	12	0.00	0.00	Manisha mufe	Active	Biomedical
	Call Closed	Biomedical/13-06-2020/4	1	Ventilator	servo's	37861	Maquet	MICU	Pre use fail	Checked and found ok	13-06-2020 22:50	13-06-2020 22:55	0.05	5	0.00	0.00	Shamila Dhage	Active	Biomedical
	Call Closed	IT/11-07-2020/1	1	Ventilator	servo's	37861	Maquet	MICU	Internet not working	Checked proxy setting.	11-07-2020 11:24	11-07-2020 11:30	0.06	6	0.00	0.00		De Active	Biomedical
	Call Closed	Biomedical/06-08-2020/2	1	Ventilator	servo's	37861	Maquet	MICU	Pre use test failed	Replaced nozzle.	06-08-2020 15:11	07-08-2020 13:00	21.49	4	4500.00	0.00	Ruchita Patil	Active	Biomedical
	Call Closed	Biomedical/08-10-2020/3	1	Ventilator	servo's	37861	Maquet	MICU	Pre-use test fail.	Check and found ok.	08-10-2020 12:20	08-10-2020 14:00	1.40	10	0.00	0.00	Biomedical Admin	Active	Biomedical
	Call Closed	Biomedical/13-10-2020/1	1	Ventilator	servo's	37861	Maquet	MICU	system volume too small.	Replaced both module nozzle	13-10-2020 16:25	13-10-2020 16:50	0.25	5	0.00	0.00	Manisha mufe	Active	Biomedical
	Call Closed	Biomedical/04-01-2021/2	1	Ventilator	servo's	37861	Maquet	MICU	Not action	Check and found	04-01-2021 09:10	04-01-2021 09:40	0.30	10	0.00	0.00	Shamila Dhage	Active	Biomedical

Equipment history-3

Spare-Part Used Details

Sr.No.	Location	Job Number	Partid	Part Name	Quantity	Price Per Item	Total	Used Date
1	Biomedical	Biomedical19-01-2021/12		Maquet ventilator Battery	1	24049.00	24049.00	20-01-2021
2	Biomedical	Biomedical04-01-2021/7		Hose Clip	1	0.00	0.00	17-02-2021
3	Biomedical	Biomedical04-01-2021/7		3 pin top (6A)	1	38.00	38.00	17-02-2021










AMC/CMC Details

Equipment	Model	Serial No	Manufactural	AMC/CMC	Cost of AMC/CMC	Last AMC/CMC	Next Due Date	Added By	IsEnable	Actions
Ventilator	servo's	37861	Maquet	AMC	0.00			Biomedical Admin	Active	  

Calibration Details


Equipment Name	Model	Serial No	Manufactural	Calibration DateTime	Calibration Details	Company Engineer Name	BioMedical Engineer	Next DueDate	Actions
----------------	-------	-----------	--------------	----------------------	---------------------	-----------------------	---------------------	--------------	---------

PM Details


Equipment Name	Model	Serial No	Manufactural	Maintenace DateTime	Company Engineer Name	Biomedical Engineer Name	Next PMDueDate	status	Action
Ventilator	servo's	37861	Maquet	5/27/2020 12:00:00 AM	BME	Biomedical Admin	27-11-2020	Active	  
Ventilator	servo's	37861	Maquet	5/27/2020 12:00:00 AM	BME Dept.	Biomedical Admin	01-11-2020	Active	  
Ventilator	servo's	37861	Maquet	11/3/2020 12:00:00 AM	Sharmila	Biomedical Admin	03-05-2021	Active	  

PM graph

Preventive Maintenance Graph



Service call and Stock Dashboard



Convenient mobile app.



Scan QR code



OUR UNIQUE FEATURES

Customizable Central Dashboard with **daily Alerts for PM, Calibration, AMC, Warranty due dates}**
Mobile App

Ai based analysis and reporting

Auto PM scheduling

Training module

E-approval's and one click Report generation
(100 % paperless)

Multiple department management.
(Biomedical, Facility & IT)

Maintenance Cost analysis

OUR UNIQUE FEATURES

Annual Budget allocation Calculations

Equipment transfer facility

Equipment utilization (ROI)

License record management

Auto calibration scheduling

Excel to software data upload

Depreciation cost calculation

SUMMARY

- 1 Complete Asset Management package
- 2 Service Call / Breakdown management
- 3 Equipment History (Total Expenses /total PM done /Total calibration)
- 4 Calibration scheduling
- 5 AMC/CMC/PM/Calibration alerts
- 6 Automatic Calculation of up-time & down-time of equipment
- 7 Gives response time of breakdown to management
- 8 Training management
- 9 Calculate cost % of total expenses / Cost analysis
- 10 Graph/pie chart representation
- 11 Location wise asset tracking
- 12 Asset transfer location to location
- 13 Inventory management location-wise
- 14 Preventive Maintenance auto scheduling
- 15 ROI Implementation.(Return On Investment).

SUMMARY

- 16 Calculates response time of complaints.
- 17 Total expenses on individual equipment.
- 18 Parts Approval (Live tracking).
- 19 Critical equipment SMS notifications.
- 20 Part wise cost.
- 21 Customize reports.
- 22 Differentiate pending service call.
- 23. Spare parts coding equipment wise.
- 24 Service call assigning (Auto or Manual).
- 25 Engineer working hours.
- 26 Customizable reports
- 27 PM/Calibration Scheduling.
- 28 Equipment transfer.
- 29 Cost analysis
- 30 Equipment analysis

SUMMARY

- 31 Stock alerts
 - 32 Equipment Cost exceeding Alerts
 - 33 Breakdown management with Colour Flag
 - 34 MIS Report(To Date - From Date)(Monthly Information system report)
 - 35 Used Part Report (Parts / Cost / Quantity)
 - 36 User Rights Managements
 - 37 MTTR Report(Mean Time To Resolution)
 - 38 Critical Equipment SMS Alerts
 - 39 Andriod application to maintain service calls
 - 40 Auto sms/email alerts to vendors
 - 41 QR code scan for complete asset details
-

Golden Bird Trading Est.



Hasan Al Majed

Managing Director

Tel.: 771 90 001 - Fax: 136 22 175 - Mobile : 39 444 791

Kingdom of Bahrain

E-mail: bbahraini@gmail.com
